



IDEXX

SNAP[®] BVDV Ag Test

Improve performance
in just 20 minutes



Test all new introductions

Backgrounder • Stocker • Feedyard

Vaccination is not enough to effectively manage BVDV

IDEXX recommends a simple three-step plan:

1. Determine PI status of all incoming cattle.
2. Remove and quarantine all positive animals.
3. Work with your veterinarian or consultant to design a vaccination program.

The SNAP® BVDV Ag Test helps identify BVDV-infected cattle in only 20 minutes to prevent the spread of BVDV

It's quick and easy to run on farm:

- Large ear-notch samples are easy to collect.
- The test is simple to run.
- Accurate ELISA results are available in 20 minutes.

Using the SNAP® BVDV Ag Test

Before you begin

Allow the test to come to room temperature (65–85°F).

Test procedure for large ear notches (≥1 cm)

1




- Dispense 1.5 mL of conjugate into a tube containing a single ear notch.
- Close the tube and swirl gently to coat ear notch.
- Let sample soak at room temperature (65–85°F) for at least 10 minutes (range: 10–60 minutes).

2



Pour liquid contents of the vial into the sample well of the SNAP® device.

3



activation circle


When sample **first** appears in the activation circle (may take 15–60 seconds), press firmly to activate. You will hear a “snap” sound. Note: Some liquid may remain in the sample well.

4



Wait 10 minutes, then read results (time range to wait: 10–15 minutes).

Result interpretation



Sample is POSITIVE if sample spot is blue. To confirm persistently infected status, run a second test 3 weeks later, using a fresh sample.



Sample is NEGATIVE if sample spot has no color development.

Test is valid if:

- Positive control spot is blue.
- Negative control spot has no color development.



Removing BVDV PI-infected cattle can reduce performance losses by as much as 42 dollars per head¹

Remove PI animals, which continually shed virus and threaten the health of the entire pen.

Improve animal health with the test that has been shown to reduce costly antibiotic use by up to 40%.²

Eliminate a major cause of a respiratory wreck in feeder cattle.³



Test with Confidence™

1-800-548-9997

idexx.com/snapbvdv



IDEXX Laboratories, Inc.

One IDEXX Drive

Westbrook, Maine 04092

United States

References

1. Hessman BE, Fulton RW, Sjeklocha DB, Murphy TA, Ridpath JF, Payton ME. Evaluation of economic effects and the health and performance of the general cattle population after exposure to cattle persistently infected with bovine viral diarrhoea virus in a starter feedlot. *Am J Vet Res.* 2009;70(1):73–85. doi:10.2460/ajvr.70.1.73
2. Sjeklocha D. Antibiotic stewardship—our approach. Presented at: 2015 NIAA Antibiotic Symposium—Stewardship: From Metrics to Management; November 4, 2015; Atlanta, GA. www.animalagriculture.org/resources/Documents/Conf%20-%20Symp/Symposiums/2015%20Antibiotics/Speaker%20Presentations/Sjeklocha,%20Dave.pdf.
3. Reggiardo C. Role of BVD in shipping fever of feedlot cattle: Case studies and diagnostic considerations. In: *Proceedings of Annual Meeting—American Association of Veterinary Laboratory Diagnosticians*. Vol 22. 1979:315–320.

© 2020 IDEXX Laboratories, Inc. All rights reserved. • 09-2497834-00
All ®/TM marks are owned by IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. The IDEXX Privacy Policy is available at idexx.com.