
## Positioning a patient for accurate vertebral heart score measurement

## Positioning is critical for accurate visualization of the heart

Guidelines:

- Place the patient in **right lateral** recumbency with no axial rotation
- Extend the forelimbs forward.
- Place the head in a neutral position.
- Take the radiograph as close to **peak inspiration** as possible.

Use weighted sandbags, foam wedges, and other positioning aids, as needed, to ensure correct positioning.




## Correct (right lateral), VHS:10.31

TIP: Find the VHS tool under IDEXX Web PACS™ Advanced Viewer.


## Why patient positioning matters

- Improper positioning can make the heart appear larger or smaller, depending on degree of rotation, and can cause cardiac chambers to appear abnormal and exaggerated in size.
- Failure to extend the forelimbs can cause poor visualization of the cranial cardiac border and cranial lung fields.
- Expiratory radiographs can mimic a diffuse pulmonary interstitial pattern.
- Improper positioning can generate misleading VHS values.



Incorrect (rotated) VHS:

10.08



Incorrect (left lateral) VHS: 9.72

Correct positioning helps to ensure an accurate clinical picture and meaningful VHS values. Visit the IDEXX Learning Center for additional training modules on radiographic positioning.

